

Instructions for Homemade Cotton Face Masks *

Supplies

- ✓ Tight-weave cotton fabric (i.e. quilting cotton). Fabric should be newly purchased (within last 12 months), washed and dried without fragrance or dyes prior to sewing. Two contrasting fabrics are needed so that the front and back are distinguishable.
- ✓ One 6" piece of 18 – 20-gauge malleable stainless steel, galvanized steel, copper or brass wire that can be formed to the nose. (not aluminum as it is too soft). **Important: Metal must be rust resistant due to sterilization process.**
- ✓ Cutting mat, ruler, rotary cutter if available
- ✓ Sewing machine, cotton thread, scissors, sewing clips/paperclips, seam ripper
(NOTE: please use pins sparingly to avoid holes)
- ✓ Jewelry pliers or needle nose pliers and wire cutters

Directions

1. Prepare the wire. Make a small loop on each end of the wire and press flat. The loops should be facing the same way. Make sure there are no rough edges that could poke the wearer or damage the mask.



2. Cut one 7" x 8" (mask body) from each contrasting tight weave cotton fabric and two pieces to 2" x 36 1/2" (tie/binding fabric). The 7" sides will be the width across the top and the bottom of the mask. The 8-inch sides will be pleated. Three pleats will make the final mask 4" deep.
3. Layer the 2 mask body pieces of contrasting fabric right sides together. Sew a 1/4" seam on the 7" sides.
4. Turn the mask right side out. Finger press and clip seams; topstitch 1/4" from both edges on the sides you just clipped.
5. Insert the prepared wire into one of the 1/4" seam allowances. Center the wire. This is now the TOP of the mask.
6. Mark 2" down from the TOP edge of the mask on both sides. Make three 1/2" accordion-style pleats along the 8" side as follows:

Make the 3 half-inch pleats starting at the 2" mark and clip each pleat. Space the pleats one after the other and clip them to make the mask 4" wide. Adjust the pleats accordingly.

The pleats do not have to be perfect as long as the mask is 4 inches deep. This is VERY important because if the mask is too short, it won't go under the chin in some cases, and if the mask is too long, it could leave too much of a gap on the side of the face.


Do not change the number of pleats as the bottom pleat is especially important to the curvature of the mask under the chin.



7. Sew the pleats down using a 1/4" seam. Repeat on the other side. When finished, the pleats will all be going down.

8. Tie/Binding

- a) Cut a 2" x 36" tight weave cotton (quilting cotton) fabric. Fold in short ends over by 1/4".
- b) Fold the fabric down middle and press to mark center.
- c) Open and fold raw edges on both sides to center.
- d) press and clip to hold it together along the length of the tie. This will eliminate the raw edges on the binding and no need for the yarn.
- e) Find the center of the tie and center of the pleat sides of mask. With centers together, insert mask into binding and clip into place.



9. Using a 3/8" seam, sew the length of the tie/binding, end to end making sure to catch the face mask. To reinforce, sew a second line of stitching on both sides of the mask where the pleats are.



FINISHED MASK



(The mask in this image is not made with cotton fabric)

The surgical mask information is provided as research information only and has not been tested for commercial use. The mask information is provided with permission for the recipient to freely use, copy and modify without restriction, subject to an obligation on the recipient to recognize University of Florida Health (UF Health) as the source of the information and design depicted in this material. The design is bare bones and materials should be locally available at hospitals that sterilize equipment. This design is for an adult mask for anyone at a higher risk for coronavirus, exposure or concern thereof. The mask information is experimental in nature and the safety or efficacy for use in humans has not been proven. It has been tested using the standard N95 fit test and can fit test pass most people. It should not be used without the wearer N95 fit testing before use. If reused it is advised to use one of the standard CDC N95 mask reuse protocols. The design, and masks built in accordance with the design, have not been approved by the FDA or NIOSH. **DISCLAIMER: THE MASK INFORMATION IS PROVIDED "AS-IS, WHERE-IS," WITHOUT REPRESENTATIONS, CONDITIONS OR WARRANTIES OF ANY KIND, WHETHER EXPRESS OR IMPLIED, INCLUDING, BUT NOT LIMITED TO, WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE, OR THAT THE USE OF THE MASK INFORMATION WILL NOT INFRINGE ANY PATENT, COPYRIGHT, TRADEMARK OR OTHER PROPRIETARY RIGHTS. THE RECIPIENT IS SOLELY RESPONSIBLE FOR DETERMINING THE APPROPRIATENESS OF USING, REPLICATING OR REDISTRIBUTING THE MASK INFORMATION AND DESIGN. IN THIS REGARD, THE RECIPIENT ASSUMES ALL LIABILITY FOR DAMAGES, OF WHATEVER NATURE AND DESCRIPTION, WHETHER IN CONTRACT OR IN TORT, WHICH MAY ARISE FROM THE USE OF THE MASK INFORMATION AND DESIGN. THE UNIVERSITY OF FLORIDA, INCLUDING ITS EMPLOYEES AND AGENTS, WILL NOT BE LIABLE TO THE RECIPIENT OR TO ANY THIRD PARTY FOR ANY LOSS, CLAIM OR DEMAND MADE BY THE RECIPIENT, OR ANY LOSS, CLAIM, DEMAND OR JUDGMENT AGAINST THE RECIPIENT BY ANY OTHER PARTY, DUE TO OR ARISING FROM THE USE OF THE MASK INFORMATION AND DESIGN BY THE RECIPIENT.**

TMH Cotton Face Mask. Version 2.4-3-20