

Cart Readiness Checklist

	The cart is plugged in and charged The cart device is powered on (Power is on the upper left corner of the Surface Pro) The cart is logged into Window correctly (no prompt for username and password) The cart is connected to the network (Look for this icon VidyoConnect is open and ready to go (See below) VidyoConnect VidyoConnect
	 ✓ simplyeconnect.health4.vidyoconnect.com Username Password Sign in
	Keep me signed in
	Username and Password are available, if "Keep me signed in" isn't checked Room Link has been given to remote party When you enter the room, make sure your speakers are on. Also, ensure your camera, and microphone are available. Look for these icons.
If a	ny of these items are red with a line through them, they are on mute.
	When you're done with your patient, please make sure to plug in the cart for the next user.